

# RevSpin RS-102, 150, & 200 Microcentrifuges



## Revolutionary Science

*Manufacturer of Precision Laboratory Equipment*

The RevSpin Microcentrifuges are stronger, faster, and contain more than twice the capacity of competing palm-sized microcentrifuges. An oversized rubber base increases the reliability of the RevSpin models by absorbing vibration and limiting bench walking. Optional PCR rotors are available for all RevSpin units. Two year warranty included.

### RS-102 specification

Capacity:	16 position rotor
Tube Size:	1.5 or 2.0mL
Speed Selections:	Off/Pulse/6,000RPM/10,000RPM
RPM:	6,000 (low), 10,000 (high)
RCF:	1,735/2,245 (low), 4,819/6,238 (high)
Power Requirements:	110-120 VAC, 60Hz, 1.5 A, grounded cord
Display:	Easy to use turn dial
Durability:	8mm oversized shaft
Safety:	Power automatically shuts if when lid is ajar



### Revolutionary Science

17319 Lake Blvd.  
Shafer, MN 55074  
Sales: (651) 674-3109  
Fax: (775) 218-7015

E-mail: [Sales@revsci.com](mailto:Sales@revsci.com)  
Website: [www.revsci.com](http://www.revsci.com)

**ISO 9001**  
REGISTERED FIRM

# RevSpin RS-102, 150, & 200

## Microcentrifuges



### Revolutionary Science

*Manufacturer of Precision Laboratory Equipment*

#### RS-150 specifications

Capacity:	16 position rotor
Tube Size:	1.5 or 2.0mL
Speed Selections:	Off/Pulse/6,000RPM/10,000RPM
RPM:	6,000 (low), 10,000 (high)
RCF:	1,735/2,245 (low), 4,819/6,238 (high)
Power Requirements:	110-230 VAC, 50-60 Hz, 1.5 A, grounded cord
Simplicity:	Two button operation
Durability:	8mm oversized shaft
Safety:	Power automatically shuts if when lid is ajar



#### RS-200 specifications

Capacity:	16 position rotor
Tube Size:	1.5 or 2.0mL
Speed Selections:	Multiple speeds up to 10,000RPM
RPM:	10,000 max
RCF:	6,238 max
Power Requirements:	110-230 VAC, 50-60 Hz, 1.5 A, grounded cord
Timer:	60 minute digital timer
Simplicity:	Easy to read digital display
Durability:	8mm oversized shaft
Safety:	Power automatically shuts if when lid is ajar



#### Revolutionary Science

13229 Saint Croix Ave.  
Lindstrom, MN. 55045  
Sales: (651) 674-3109  
Fax: (775) 218-7015

E-mail: [Sales@revsci.com](mailto:Sales@revsci.com)  
Website: [www.revsci.com](http://www.revsci.com)

**ISO 9001**  
REGISTERED FIRM