

## **ELECTROPHORESIS POWER SUPPLY**

**\$281.00**



Electrophoresis Power Supply works with the [RevBox](#) for running routine electrophoresis gels. It is capable of running two gels simultaneously at the same or different voltages.



---

**SKU:** RS-PS-75

**Category:** [Electrophoresis](#)

## **PRODUCT DESCRIPTION**

### **ELECTROPHORESIS POWER SUPPLY**

#### **Specifications**

Max Power: 52.5 watts

Max Current: 350mA

Input Voltage: 100-240 VAC

Output Voltages: 0 volts, 15 volts, 75 volts, 150 volts

Output Jack Sets: Run 2 gel boxes simultaneously

Fuse: Glass encapsulated

Net Weight: 0.68kg (1.5lbs)

Dimensions: 7cm tall x 12.7cm wide x 10.2cm deep

2.8" tall x 5" wide x 4" deep

#### **Features**

Glass encapsulated fuse

Two output jacks

Highly efficient design provides accurate voltage output.

Red warning LED lights up when output current exceeds safe level

Easy to use slide switches

Detachable grounded cord

2 year warranty

